

Standard Operating Procedure (SOP)	
SOP Title	Laboratory protocol for testing and rebalancing for the new primer panels.
SOP number	GbS02a – 96 well
SOP Version	2

## 1. SCOPE

This SOP describes the procedure for determining the optimum concentration of each primer pair in a genotyping by sequencing (GbS) reaction for even amplification across the individual targets within a panel.

Before batch processing of samples, a new GbS design requires test amplification of DNA samples using an equimolar pool of primers. This is because the performance of individual primer pairs within a multiplex cannot be predicted *in silico* so must be determined empirically. The test libraries are then Ampure SPRI size selected and purified and following quality control checks pooled before sequencing as a 150PE MiSeq run in an identical manner to that of a production GbS run.

The resultant read data are analysed to provide a result of the level of coverage of each target as a percentage of all on target reads.

Concentrations of each primer pair are subsequently calculated and adjusted to compensate for over- and underperforming primers. A rebalanced pool is prepared and aliquoted for future use as a working panel.

## 2. MATERIALS REQUIRED:

2.1 REAGENTS AND MATERIALS		
<b>Primers</b>		
Stock oligonucleotide GbS panel @ 500µM - GCR1, GCR2, & Speciation (SPEC)		
i7 & i5 lyophilised GbS primers		
<b>Reagents</b>	<b>Supplier</b>	<b>Cat. No.</b>
QIAGEN Multiplex PCR polymerase (x2)	Qiagen	206145
Ampure XP beads (SPRI) - Beckman Coulter	Fisher Scientific	NC9959336
Nuclease free water (NFW) - Ambion	Fisher Scientific	AM9930
Ethanol (75%)		
T0.1E buffer (10 mM Tris, 0.1 mM EDTA, pH 8.0)		
Buffer Elution Buffer (Buffer EB, 10 mM Tris-Cl, pH 8.5)	Qiagen	19086
Genomic DNA High Sensitivity D1000 ScreenTape	Agilent	5067-5584
Genomic DNA High-sensitivity reagents	Agilent	5067-5585
KAPA SYBR FAST ABI Prism qPCR Kit, KAPA SYBR FAST qPCR Master Mix (2X) Library Quantification DNA Standards 1 – 6 Library Quantification Primer Premix (10X)	Sigma-Aldrich	KK4605
<b>Samples</b>	<b>Supplier</b>	<b>Cat. No.</b>
Test DNA sample. 2µL sWGA/PEP amplified DNA		
<b>Materials</b>	<b>Supplier</b>	<b>Cat. No.</b>
PCR machine		
Microplate heat sealer (PX1)	Bio-Rad	#1814000
Foil pierceable plate seal (heat sealer)	Bio-Rad	#1814040
Tissue wipes (Azowipes)		
Lo-bind Eppendorf tubes (1.5mL)		
Ice tray/cooling block		
96 well qPCR plate, skirted (AB2800)	Fisher Scientific	10032013
Adhesive PCR Film Polyester plate seals	Thermo Fisher Scientific	AB0558
Tube/ Microplate vortexer		
Magnetic tube rack (Invitrogen DynaMag)	Thermo Fisher Scientific	12321D
Non-Magnetic tube rack		
Centrifuge		
Agilent TapeStation System	Agilent	
Agilent – tube vortexer (IKA MS3)	included with TapeStation system	
Agilent consumables (loading tips, optical tube, optics cap, optical 96 well plate and foil seals)	Specific to TapeStation system	
qPCR Instrument (ROCHE LightCycler 480II)	ROCHE	
White 96 well qPCR plate	StarLab	I1402-9909
Optically clear plate seal	StarLab	E2796-9795
8-channel Adjustable Multichannel (0.5-10µL)		
12-channel Adjustable Multichannel (10-100µL)		
Pipettes (1000, 200, 20 and 1 µL)		
Filtered tips of various volumes (compatible with pipettes and multichannel)		

## 2.2 PCR\_1 PRIMER PANEL

GbS primers (100nmol scale synthesis from IDT, supplied re-suspended in T0.1E at 500 $\mu$ M in a 96 deep well plates); it is assumed that the primers have been ordered with paired forward and reverse primers in separate plates (e.g. forward primers in Plate 1 A1, B1..., reverse primers in Plate 2 A1, B1...)

PCR is sensitive to inhibitors, which may be introduced due to a DNA source or introduced via DNA extraction. A review article can be found here for further details.

Oyola, S.O., Ariani, C.V., Hamilton, W.L. et al. Whole genome sequencing of Plasmodium falciparum from dried blood spots using selective whole genome amplification. Malar J 15, 597 (2016).

## 2.3 PCR\_2 Primer Panel

Pre-made PCR\_2 tag plates as produced as in Tag plate preparation SOP.

## 2.4 PRIMER DILUENT

The primer diluent (T1E0.1 +0.01% v/v Triton X-100) is used in section 9 for the diluting of the pooled primers to 40nM before aliquoting.

To create the primer diluent do the following: Dilute 10mM Tris/1mM EDTA pH 8 (T10E1) 1in10 with 10mM Tris pH 8 (e.g. Qiagen EB) then spiking in the 0.01%v/v Triton X-100 detergent.

- Add 5mL of 10x TE (pH 8) to a 50mL falcon tube.
- Then add 45mL of EB (10mM Tris).
- Dispense 5 $\mu$ L of neat Triton X-100 onto the inner surface of the tube slowly, using a P10 pipette, i.e. not into the diluent.
- Vortex thoroughly until all the Triton X-100 is in solution.

## 2.5 SAMPLES AND BATCH SIZE

Please note that sample selection is very important. For the *P. falciparum* primer panel only *P. falciparum* positive samples must be used. Furthermore, for the *P. vivax* primer panel only *P. vivax* positive samples must be used.

Input DNA does not require prior quantification or normalisation. Normalisation of library yields are achieved during the PCR itself.

Each plate must contain at least one negative control (i.e. **2µl of T0.1E**).

Test DNA samples (of a representative type/quality/concentration to that expected to be supplied in production samples):

- All samples must have a volume of 2µL.
- 95 of the 96 wells must have sample (96<sup>th</sup> well is the negative control).
- Either a minimum of 95 individual samples can be used or 24 samples in multiples across the plate.

A large number of negative/empty wells will generate significant amounts of smaller PCR artefacts (primer dimers), which are problematic to fully remove by size selection and will preferentially sequence.

Therefore, the plate must have all 96 wells used for balancing.

### 3. LIBRARY PREPARATION

Primers
Plates containing GRC1, GRC2 and Speciation GbS primers. Forward and Reverse primers for each set.

PCR Reagents
QIAGEN Multiplex PCR polymerase (x2)
Stock oligonucleotide GbS panel @ 500µM - GCR1, GCR2, & Speciation (SPEC)
Nuclease free water (NFW)
T0.1E buffer

PCR Equipment		
PCR Machine	Foil pierceable plate seal	1000 µl pipette
96 well qPCR plate	Centrifuge	200 µl pipette
Ice tray/cooling block	Tissue wipes (Azowipes)	20 µl pipette
Microplate heat sealer		

3.1 PRIMER POOLING	
See Appendix 4.	
3.1.01	Clear out a sufficient work area and wipe down the lab bench surface and pipettes with tissue wipes (Azowipes) or similar.
3.1.02	Defrost the forward and reverse primer plates for each target (GRC1, GRC2 & Speciation @ 500µM).
3.1.03	Once defrosted, gently vortex the plate to ensure the primers are fully re-suspended.
3.1.04	Centrifuge the plate to collect all the liquid at the bottom of the wells.
3.1.05	Add 7.5µL from each forward primer well into the equivalent wells of a new 96-well PCR plate (equivolume primer plate).
3.1.05	Add 7.5µL from each reverse primer well into the equivalent wells of the “equivolume primer plate” The resultant wells of the “equivolume primer plate” should contain 15µL of a 1:1 mixture of forward and reverse primers (at 250µM).
3.1.06	The “equivolume primer plate” can now be used as the stock for both the “equivolume” and “rebalanced” primer pools in order to minimise accumulated pipetting errors.

Equivolume pool	
3.1.07	Centrifuge the plate of paired primers.
3.1.08	Using a 12-channel remove 2µL from each primer pair well and pool into a single column of empty wells on a new 96 well plate.
3.1.09	Then using a single channel combine the subpools together to produce the “equivolume pool” into an appropriately labelled Eppendorf tube.
3.1.10	<p>Calculate the concentration of the individual primers in the pool. Label the Eppendorf with calculated concentration.</p> <p>Concentration of each individual primers = <math>250\mu\text{M} \div (\text{no. primer pairs})</math></p> <p>GRC1 = <math>250/68 = 3.67\mu\text{M}</math></p> <p>GRC2 = <math>250/67 = 3.73\mu\text{M}</math></p> <p>Speciation = <math>250/2 = 125\mu\text{M}</math></p>
Working pool	
3.1.11	<p><b>Prepare a “working pool” aliquot of the stock at the working concentration (40nM) in a Lo-bind Eppendorf (dilute in T0.1E).</b></p> <p>Dilution factor = <math>([\text{stock}]/40\text{nM})</math></p> <p>GRC1 = <math>(3.67\mu\text{M} \times 1000)/40\text{nM} = 91.75</math></p> <p>Pooled primer pool = <math>400/91.75 = 4.4\mu\text{L}</math></p> <p>Diluent (T0.1E) = <math>400 - 4.4 = 395.6\mu\text{L}</math></p> <p>GRC2 = <math>(3.73\mu\text{M} \times 1000)/40\text{nM} = 93.25</math></p> <p>Pooled primer pool = <math>400/93.25 = 4.3\mu\text{L}</math></p> <p>Diluent (T0.1E) = <math>400 - 4.3 = 395.7\mu\text{L}</math></p> <p>Speciation (pre-dilution) = <math>1:25 = 1\mu\text{L} + 24\mu\text{L Diluent (T0.1E)}</math></p> <p>= <math>(125/25) = 5\mu\text{M}</math></p> <p>Speciation = <math>(5\mu\text{M} \times 1000)/40\text{nM} = 125</math></p> <p>Pooled primer pool = <math>400/125 = 3.2\mu\text{L}</math></p> <p>Diluent (T0.1E) = <math>400 - 3.2 = 396.8\mu\text{L}</math></p>
3.1.12	Aliquot the “equivolume pool” and “working pool” into 20µL aliquots and store at -20°C. After this proceed directly to PCR_1 (section 3.2).
<p><b>Plate Storage:</b> Seal the Forward primer, Reverse primer and mixed “equivolume primer plate” using the plate sealer heat seal the lid to the plate (175 °C for 5 seconds).</p> <p>NOTE: It is important to use these exact settings, as too high a temperature will melt the top of the wells and too low will lead to incomplete sealing and evaporation during PCR.</p> <p style="text-align: center;"><b>Plates can be stored at -20°C until required.</b></p>	

### 3.2 METHODOLOGY PCR\_1

During PCR\_1, genomic priming takes place via a low number of PCR cycles utilising a long annealing stage, which allows the highly dilute primers to specifically anneal to their genomic target (see Appendix 5).

#### PCR Reagents

Primer Panel “equivolume working pool” (40nM)

QIAGEN Multiplex PCR polymerase (x2)

Buffer EB

Nuclease free water (NFW)

#### Samples

Test DNA samples (sWGA or PEP DNA - 2µL)

#### PCR Equipment

PCR machine	Foil pierceable plate seal	1000 µl pipette
96 well qPCR plate	Centrifuge	200 µl pipette
Ice tray/cooling block	Tissue wipes (Azowipes)	20 µl pipette
Microplate heat sealer		1µL pipette

#### Methodology PCR\_1 – Pre-PCR room

3.2.01	Clear out a sufficient work area and wipe down the lab bench surface and pipettes with tissue wipes (Azowipes) or similar.
3.2.02	Prepare ice tray/cooling block for use during PCR_1.
3.2.03	Aliquot <b>2µL</b> of <b>each</b> of the <b>test DNA samples</b> – a minimum 24 individual samples in multiples spaced over the 96 well plate (please ensure you meet minimum plate coverage - see section 2.5)
3.2.04	Add 2µL of buffer EB to the one or two empty wells (negative control wells).
3.2.05	Spin the plate down in the centrifuge to ensure that all liquid is collected at the bottom of the wells.

**Note** - Start PCR\_1 setup first thing in the morning to enable processing through to PCR\_2 within a working day.

3.2.06	<p>Make the PCR_1 mastermix for the number of plates being processed as shown below (Table 1) in a Lo-bind Eppendorf tubes (1.5mL). The “x110” will be sufficient for a single 96 well plate.</p> <p>N.B. It is important to pipette up and down when adding the primer panel to ensure the entire volume is transferred into the master mix.</p> <table border="1" data-bbox="354 430 1279 712"> <thead> <tr> <th></th> <th></th> <th>1 plate</th> <th>2 plates</th> </tr> <tr> <th>Components:</th> <th>Per well (µL)</th> <th>x110 (µL)</th> <th>X250 (µL)</th> </tr> </thead> <tbody> <tr> <td>NFW</td> <td>3.433</td> <td>377.6</td> <td>855</td> </tr> <tr> <td>Qiagen MM (x2)</td> <td>5.500</td> <td>605</td> <td>1375</td> </tr> <tr> <td>“working pool” (40nM) 1/133</td> <td>0.068</td> <td>7.48</td> <td>17</td> </tr> <tr> <td><b>TOTAL</b></td> <td><b>9</b></td> <td><b>990</b></td> <td><b>2247</b></td> </tr> </tbody> </table> <p>Table 1: 2µL DNA in sample plate; volumes of mastermix for 1 or 2 plates including excess.</p>			1 plate	2 plates	Components:	Per well (µL)	x110 (µL)	X250 (µL)	NFW	3.433	377.6	855	Qiagen MM (x2)	5.500	605	1375	“working pool” (40nM) 1/133	0.068	7.48	17	<b>TOTAL</b>	<b>9</b>	<b>990</b>	<b>2247</b>
		1 plate	2 plates																						
Components:	Per well (µL)	x110 (µL)	X250 (µL)																						
NFW	3.433	377.6	855																						
Qiagen MM (x2)	5.500	605	1375																						
“working pool” (40nM) 1/133	0.068	7.48	17																						
<b>TOTAL</b>	<b>9</b>	<b>990</b>	<b>2247</b>																						
3.2.07	Gently mix the master mix by vortexing and pulse spin down.																								
3.2.08	Add <b>9µL</b> of PCR_1 master mix to each well of a plate and mix.																								
3.2.09	Apply a foil pierceable plate seal on top of the plate, red strip facing upwards.																								
3.2.10	Using the plate sealer heat seal the lid to the plate (175 °C for 5 seconds). NOTE: It is important to use these exact settings, as too high a temperature will melt the top of the wells and too low will lead to incomplete sealing and evaporation during PCR.																								
3.2.11	Centrifuge at 1000g for 20 seconds at 4°C																								
3.2.12	<p>Transfer to a PCR machine, tighten the lid and start the required program using the parameters below (Table 2: PCR_1):</p> <table border="1" data-bbox="453 1361 1270 1617"> <thead> <tr> <th colspan="4">Standard PCR_1 – GRC1, GRC2, SPEC</th> </tr> </thead> <tbody> <tr> <td>Step 1</td> <td>95°C</td> <td>15:00 min</td> <td rowspan="5">} Total of 5 cycles</td> </tr> <tr> <td>Step 2</td> <td>95°C</td> <td>00:20 sec</td> </tr> <tr> <td>Step 3</td> <td>51°C</td> <td>40:00 min</td> </tr> <tr> <td>Step 4</td> <td>60°C</td> <td>03:00 min</td> </tr> <tr> <td>Step 5</td> <td>4°C</td> <td>∞</td> </tr> </tbody> </table> <p>Table 2: PCR_1 programs and cycling conditions.</p>	Standard PCR_1 – GRC1, GRC2, SPEC				Step 1	95°C	15:00 min	} Total of 5 cycles	Step 2	95°C	00:20 sec	Step 3	51°C	40:00 min	Step 4	60°C	03:00 min	Step 5	4°C	∞				
Standard PCR_1 – GRC1, GRC2, SPEC																									
Step 1	95°C	15:00 min	} Total of 5 cycles																						
Step 2	95°C	00:20 sec																							
Step 3	51°C	40:00 min																							
Step 4	60°C	03:00 min																							
Step 5	4°C	∞																							
3.2.13	Prepare ice tray/cooling block for use after PCR_1.																								
3.2.14	Once the PCR Machine reaches “4°C” the PCR_1 program is complete, <b><u>immediately place on ice/cold block.</u></b>																								
3.2.15	Proceed to next stage (see section 3); it is recommended that the plate is left at 4°C for no longer than three hours.																								



### 3.3 METHODOLOGY PCR\_2

The second PCR involves no additional reagents or clean up after the PCR\_1. A volume of PCR\_1 will be transferred to a dried down dual indexed tag plate. These indexed tag primers anneal to and extend from the 5'-tails incorporated during PCR\_1. This enables the capability to multiplex up to 1536 samples per sequencing batch (see Appendix 5).

#### PCR Reagents

PCR\_2 Tag plates (i5 and i7 lyophilised primers) – **400nM**

Nuclease free water (NFW)

#### PCR Equipment

Ice tray/cooling block	Foil pierceable plate seal	Centrifuge
PCR machine	Tube/ Microplate vortexer	Multi-channel (10µL)
Microplate heat sealer	Tissue wipes (Azowipes)	20µl pipette

#### Methodology PCR\_2 plate preparation

It is important to continue to PCR\_2 as soon as possible after PCR\_1 has finished.

3.3.01	Take the completed PCR_1 plates and PCR_2 tag plates and centrifuge at 2000g for 1 minute (at 4°C) and place back on ice tray/cooling block.
3.3.02	Place the PCR_1 plates containing the test reactions and a pre-prepared PCR_2 96 well tag plates on ice. Make a note of which i5 primers are used, this information is required at the submission stage.
3.3.03	Remove the seals from the PCR and Tag plates.
3.3.04	Using a 10µL multichannel pipette (either 8 or 12 channel), transfer 10µL of the PCR_1 reaction volumes into the equivalent well of the PCR_2 tag 96 well plates and mix thoroughly (mix by pipetting 10 times). Ensure that the dried tag primers in the PCR_2 plate are fully re-suspended.
3.3.05	Apply a foil pierceable plate seal on top of the plate, red strip facing upwards.
3.3.06	Using the plate sealer heat seal the lid to the plate (175 °C for 5 seconds). NOTE: It is important to use these exact settings, as too high a temperature will melt the top of the wells and too low will lead to incomplete sealing and evaporation during PCR.
3.3.07	Keep created PCR_2 Tag plate on ice and transfer to a centrifuge. Centrifuge plates at 2000g for 1 minute (at 4°C). Immediately place back on ice.

Methodology PCR_2 – Post-PCR room																																		
<b><u>PLEASE NOTE AT THIS STAGE IF POSSIBLE IT IS ADVISABLE TO PROCEED TO THE POST-PCR ROOM</u></b>																																		
3.3.08	Leave the plate on ice.																																	
3.3.09	<p>Pre-heat the PCR machine by tightening the lid and starting the PCR machine with the required program using the parameters below (Table 3: PCR_2).</p> <table border="1" style="margin-left: auto; margin-right: auto;"> <thead> <tr> <th colspan="5">GRC1 / GRC2 / SPEC (Subcycling PCR_2)</th> </tr> </thead> <tbody> <tr> <td>Step 1</td> <td>95°C</td> <td>∞</td> <td></td> <td></td> </tr> <tr> <td>Step 2</td> <td>95°C</td> <td>00:20 sec</td> <td></td> <td></td> </tr> <tr> <td>Step 3</td> <td>68°C</td> <td>00:15 sec</td> <td rowspan="2">} Total of 4 cycles</td> <td rowspan="2">} Total of 31 cycles</td> </tr> <tr> <td>Step 4</td> <td>60°C</td> <td>00:15 sec</td> </tr> <tr> <td>Step 5</td> <td>68°C</td> <td>03:00 min</td> <td></td> <td></td> </tr> <tr> <td>Step 6</td> <td>4°C</td> <td>∞</td> <td></td> <td></td> </tr> </tbody> </table> <p>Table 3: PCR_2 programs and cycling conditions</p>	GRC1 / GRC2 / SPEC (Subcycling PCR_2)					Step 1	95°C	∞			Step 2	95°C	00:20 sec			Step 3	68°C	00:15 sec	} Total of 4 cycles	} Total of 31 cycles	Step 4	60°C	00:15 sec	Step 5	68°C	03:00 min			Step 6	4°C	∞		
GRC1 / GRC2 / SPEC (Subcycling PCR_2)																																		
Step 1	95°C	∞																																
Step 2	95°C	00:20 sec																																
Step 3	68°C	00:15 sec	} Total of 4 cycles	} Total of 31 cycles																														
Step 4	60°C	00:15 sec																																
Step 5	68°C	03:00 min																																
Step 6	4°C	∞																																
3.3.10	Once the temperature has been reached, loosen the lid and place the plate on the PCR machine and then tighten the lid.																																	
3.3.11	Select continue/skip to move to Step 2 to begin the PCR_2 program.																																	
3.3.12	The PCR can remain on the PCR machine overnight at 4°C without any impact on the downstream steps.																																	
3.3.13	Once the PCR is complete, you can either move directly to pooling or plates may be stored at -20°C for up to one week.																																	
<p>As the Ampure XP beads size selection is extremely sensitive to errors in pipetting volumes, it is recommended to ensure that the pool volume carried forward into size selection is a constant 100µL.</p>																																		

### 3.4 POOLING, PURIFICATION & SIZE SELECTION.

Size selection takes place via a 0.75x SPRI. Under these conditions the majority of small primer dimers and unincorporated primers <170bp will remain in the supernatant whilst the desired amplicons bind onto the beads (see Appendix 6).

#### Reagents

Ampure XP beads (SPRI)

75% ethanol

Buffer EB

#### Equipment

Centrifuge	1.5ml Lo-bind Eppendorf tubes	1000µL pipette
Magnetic tube rack	Tube/ Microplate vortexer	200µL pipette
Non-Magnetic tube rack		100µL pipette
Tissue wipes (Azowipes)		20µl pipette

#### Pooling of PCR products

3.4.01	Remove Ampure XP beads from fridge and allow to warm to room temperature (approx. 30min).
3.4.02	Take plates off the PCR machine and place on ice.
3.4.03	Label 3 new 1.5ml Lo-bind Eppendorf tubes and label them <ul style="list-style-type: none"> <li>TARGET pool stock (e.g. GRC1 pool stock).</li> </ul>
3.4.04	Transfer the contents of the 96 well plate from PCR2 into a fresh 96 well plate using an multi-channel pipette (either 8 or 12). Mix and transfer to the labelled 1.5ml Lo-bind Eppendorf tube (GbS pool stock) using a pipette and discard the PCR_2 plate. <b>NOTE: When pooling keep each Target plate separate.</b>
3.4.05	At this stage, it is possible to view the size profile of the pool by analysing on an Agilent TapeStation (Agilent High Sensitivity D1000 Kit). An example size profile is shown in Section 4 & Appendix 7.

Size selection and Clean-up	
3.4.06	Vortex Ampure XP beads to re-suspend
3.4.07	Accurately transfer 100µL of pooled PCR product into a new 1.5ml Lo-bind Eppendorf. Make sure the tube is clearly labelled with target name (e.g. GRC1).
3.4.08	Then add 75µL of Ampure XP beads and close the lids.
3.4.09	Vortex to mix and pulse spin down to collect all of the liquid at the bottom of the tube.
3.4.10	Place on a non-magnetic rack at room temperature for 5 minutes.
3.4.11	Transfer tube to a magnetic tube rack until clear (approx. 3-4min).
3.4.12	Once the supernatant is clear and all beads are captured on the side of the magnet, remove and discard the supernatant. Be careful not to disturb or aspirate the beads.
3.4.13	Using a 1000µL pipette and keeping the tube on the rack, wash the beads by adding 700µL of 75% ethanol, being careful not to disturb the beads captured on the magnet.
3.4.14	After 30 seconds, carefully remove and discard the ethanol
3.4.15	Repeat the ethanol wash once more (3.4.13 and 3.4.14)
3.4.16	Close the tube lid and pulse centrifuge to draw final residual ethanol to the bottom and remove.
3.4.17	With the lid open, air dry on a magnetic rack at room temperature for approximately 2 minutes.
3.4.18	Off magnet, re-suspend the beads in 105µl of Buffer EB and vortex to mix.
3.4.19	Incubate <b>off-magnet</b> for 5 minutes to release the target product into solution. Spin down if any splashes have occurred on the side of the tube.
3.4.20	Transfer to the magnetic rack.
3.4.21	Label 3 new 1.5ml Lo-bind Eppendorf tubes (1 per Target): First TARGET size selections and the relevant I5 primer (e.g. 1 <sup>st</sup> GRC1-T5-1).
3.4.22	Once the supernatant is clear (approx. 3-4min) and beads are captured on the magnet, transfer all of the supernatant containing the purified size selected PCR products to the Lo-bind 1.5ml Eppendorf tube labelled "1st TARGET selection".

Final size selection on ELUTION (see Appendix 4)	
3.4.23	Vortex Ampure XP beads to re-suspend.
3.4.24	Transfer 100µl of the “First TARGET selections” to a new Eppendorf. Label the new Eppendorf with “Second TARGET selection” and the relevant I5 primer (e.g. 2 <sup>nd</sup> GRC1-T5-1).
3.4.25	Add 75µL of Ampure XP beads to the “2nd TARGET selection” Eppendorf and close the lid.
3.4.26	Vortex to mix and pulse spin down to collect all of the liquid at the bottom of the tube.
3.4.27	Place on a non-magnetic rack at room temperature for 5 minutes.
3.4.28	Transfer tube to a magnetic tube rack until clear (approx. 3-4min).
3.4.29	Once the supernatant is clear and all beads are captured on the side of the magnet, remove and discard the supernatant. Be careful not to disturb or aspirate the beads.
3.4.30	Using a 1000µL pipette and keeping the tube on the rack, wash the beads by adding 700µL of 75% ethanol, being careful not to disturb the beads captured on the magnet.
3.4.31	After 30 seconds, carefully remove and discard the ethanol
3.4.32	Repeat the ethanol wash once more (3.4.30 and 3.4.31).
3.4.33	Close the tube lid and pulse centrifuge to draw final residual ethanol to the bottom and remove.
3.4.34	With the lid open, air dry on a magnetic rack at room temperature for approximately 2 minutes. NOTE: At this point, the size selected DNA will be dried onto the beads.
3.4.35	Re-suspend beads with 22µL of Buffer EB by pipetting up and down.
3.4.36	Transfer to non-magnetic rack and incubate “off magnet” for 2 minutes.
3.4.37	Transfer the tube containing eluted size selected libraries to the magnetic rack.
3.4.38	Allow supernatant to clear as beads are captured on the magnet (approx. 3-4min).
3.4.39	Transfer supernatant containing the purified size selected PCR products to 1.5ml Lo-bind Eppendorf tube labelled “TARGET Pool” (e.g. GRC1 pool).
3.4.40	The <b>remainder of the pooled PCR stocks</b> should be stored at -20°C.

## 4. LIBRARY QUALITY CONTROL

### Agilent TapeStation

The efficiency of the clean-up and size selection can be assessed by running 1 $\mu$ L of the “TARGET Pool” on an Agilent TapeStation using an Agilent High Sensitivity Kit.

Note: small peak at ~60bp, is the surplus primers and the small peak at ~160bp, is the result of primer dimers between first round genomic PCR primers. Due to their size, these will preferentially cluster during sequencing and result in loss of on-target data.

If the peaks are still present post clean-up and size selection then repeat Ampure XP beads size selection until the ~60bp and ~160bp peaks are barely visible by TapeStation analysis.


See Agilent TapeStation SOP.

## 5. LIBRARY QUANTIFICATION

Library quantification is performed by amplifying the set of five pre-diluted DNA Standards and diluted library samples by qPCR, using the KAPA SYBR FAST qPCR mastermix and primers targeting the Illumina P5 and P7 flow cell oligo sequences. The average C<sub>q</sub> score for each DNA Standard is plotted against log<sub>10</sub> (concentration in pM) to generate a standard curve. The concentrations of diluted library samples are then calculated against the standard curve, using absolute quantification.

See qPCR Quantification SOP.

## 6. MISEQ POOL CALCULATOR

Pool calculator (see Appendix 8)																																																																																																															
6.0.01	<p>Open the excel spreadsheet “flexible pooling calculator”.</p>  <table border="1"> <thead> <tr> <th>Pool no.</th> <th>[Pool diluted] by KAPA qPCR/pM</th> <th>Pool dilution factor (1in x)</th> <th>[Pool corrected] /nM</th> <th>Tags used</th> <th>Panel name</th> <th>⇒ Tags × Targets</th> <th>Mean [tag-target]/pM</th> <th>X relative to required final [tag-target]</th> <th>Volume needed for 4nM pool</th> <th>Buffer EB</th> </tr> </thead> <tbody> <tr> <td>MA_GRC1_1:1k</td> <td>135.50</td> <td>1000</td> <td>174.9885714</td> <td>96</td> <td>PFA_GRC1</td> <td>6624</td> <td>26.42</td> <td>88.13</td> <td>5.47</td> <td></td> </tr> <tr> <td>MA_GRC2_1:1k</td> <td>193.00</td> <td>1000</td> <td>249.2457143</td> <td>96</td> <td>PFA_GRC2</td> <td>6528</td> <td>38.18</td> <td>127.37</td> <td>3.78</td> <td></td> </tr> <tr> <td>MA_SPEC_1:1k</td> <td>21.47</td> <td>1000</td> <td>27.72266667</td> <td>96</td> <td>PFA_SPEC</td> <td>192</td> <td>144.39</td> <td>481.68</td> <td>1.00</td> <td></td> </tr> <tr> <td>4</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>5</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td>6</td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> <td></td> </tr> <tr> <td colspan="6"></td> <td>TOTAL</td> <td>13344</td> <td>✓</td> <td>471.43</td> <td>Buffer EB</td> </tr> <tr> <td colspan="6"></td> <td>⇒ Required final [tag-target] in the pool</td> <td>0.300</td> <td>pM</td> <td></td> <td></td> </tr> <tr> <td colspan="6"></td> <td>MAX</td> <td>55000</td> <td></td> <td></td> <td></td> </tr> </tbody> </table>	Pool no.	[Pool diluted] by KAPA qPCR/pM	Pool dilution factor (1in x)	[Pool corrected] /nM	Tags used	Panel name	⇒ Tags × Targets	Mean [tag-target]/pM	X relative to required final [tag-target]	Volume needed for 4nM pool	Buffer EB	MA_GRC1_1:1k	135.50	1000	174.9885714	96	PFA_GRC1	6624	26.42	88.13	5.47		MA_GRC2_1:1k	193.00	1000	249.2457143	96	PFA_GRC2	6528	38.18	127.37	3.78		MA_SPEC_1:1k	21.47	1000	27.72266667	96	PFA_SPEC	192	144.39	481.68	1.00		4											5											6																	TOTAL	13344	✓	471.43	Buffer EB							⇒ Required final [tag-target] in the pool	0.300	pM									MAX	55000			
Pool no.	[Pool diluted] by KAPA qPCR/pM	Pool dilution factor (1in x)	[Pool corrected] /nM	Tags used	Panel name	⇒ Tags × Targets	Mean [tag-target]/pM	X relative to required final [tag-target]	Volume needed for 4nM pool	Buffer EB																																																																																																					
MA_GRC1_1:1k	135.50	1000	174.9885714	96	PFA_GRC1	6624	26.42	88.13	5.47																																																																																																						
MA_GRC2_1:1k	193.00	1000	249.2457143	96	PFA_GRC2	6528	38.18	127.37	3.78																																																																																																						
MA_SPEC_1:1k	21.47	1000	27.72266667	96	PFA_SPEC	192	144.39	481.68	1.00																																																																																																						
4																																																																																																															
5																																																																																																															
6																																																																																																															
						TOTAL	13344	✓	471.43	Buffer EB																																																																																																					
						⇒ Required final [tag-target] in the pool	0.300	pM																																																																																																							
						MAX	55000																																																																																																								
6.0.02	For each sample, enter their name under the “pool No.”																																																																																																														
6.0.03	For each sample, enter the quantity mean, as determined by the qPCR under the “[pool diluted] by KAPA qPCR/pM”.																																																																																																														
6.0.04	For each sample enter the dilution factor under the “Pool dilution factor (1in x)”.																																																																																																														
6.0.05	For each sample, enter 96 under “Tags used”.																																																																																																														
6.0.06	<p>For each sample, select the correct panel from the drop down menu in “panel name”.</p> <ul style="list-style-type: none"> <li>• PFA_GRC1</li> <li>• PFA_GRC2</li> <li>• PFA_SPEC</li> </ul>																																																																																																														
6.0.07	<p>Number of Tags and Targets is automatically calculated: Check that the “TOTAL” (e.g. 13,344) is less than the “MAX” “tags x targets” (i.e. 55,000). If the number exceeds the “MAX” (55,000), then adjust the constituents of the pool to reduce the number. If the number is lower than the “MAX” (55,000) then a “✓” will appear and proceed to preparation of the 4nM pool.</p>																																																																																																														
6.0.08	Completing these fields will give you the total volume needed for your pool (blue arrow in figure above see 6.0.01).”																																																																																																														
6.0.09	Label a new Eppendorf with “GbS MiSeq Pool” and add the quantity of each pool and Buffer EB as calculated by the excel spreadsheet.																																																																																																														

## 7. MISEQ SUBMISSION

MiSeq Submission	
7.0.01	Material ready for sequencing. See Illumina MiSeq Sequencing SOP.

## 8. DATA ANALYSIS

Data analysis	
8.0.01	Following the MiSeq run proceed to testing and balancing data analysis SOP.



## APPENDIX 1 – GRC1 PRIMERS

Well	Name	Forward sequence	Reverse Sequence
A1	CRT_220	ACACTCTTTCCCTACACGACGCTCTCCGATCTATCTTTGAAACACAAGAAGmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTATTTCCCTTGTCATGTTTGAmAA
B1	CRT_326	ACACTCTTTCCCTACACGACGCTCTCCGATCTGAGCATGGGTAAAGAAGCTTAmUA	TCGGCATTCTGCTGAACCGCTCTCCGATCTCCTCTGTATGTATCAACGTmUT
C1	CRT_371	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGGTACAACGTATCATATTTmUA	TCGGCATTCTGCTGAACCGCTCTCCGATCTACGAACAAGCCATTTGATATmUA
D1	CRT_72_74_75_76	ACACTCTTTCCCTACACGACGCTCTCCGATCTTAAACAGATGGCTCACGTTmUA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGAGTTTCGGATGTTACAAAAmCT
E1	DHFR_16_51_59	ACACTCTTTCCCTACACGACGCTCTCCGATCTGTTTTCGATATTTATGCCATAmUG	TCGGCATTCTGCTGAACCGCTCTCCGATCTACATTATCCACAGTTTCTTTmGT
F1	DHFR_306	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGATGATGATGAAGAAGAAGAmUG	TCGGCATTCTGCTGAACCGCTCTCCGATCTTCGATCACTTTGTTTATTTmCA
G1	DHPS_436_437	ACACTCTTTCCCTACACGACGCTCTCCGATCTTTTTGTTGAACCTAAACGTmGC	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATAAATTGGTTTCGCATCmCA
H1	K13_resistance_1	ACACTCTTTCCCTACACGACGCTCTCCGATCTATGAATTTAGAATTCGCCAmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTCCATATGCCTTATTAGAAGmCT
A2	K13_resistance_3	ACACTCTTTCCCTACACGACGCTCTCCGATCTCATAGCTGATGATCTAGGmGG	TCGGCATTCTGCTGAACCGCTCTCCGATCTGAGGTGTATGATCGTTTAmAG
B2	K13_resistance_5	ACACTCTTTCCCTACACGACGCTCTCCGATCTTAGACATAGGTGTACACATAmCG	TCGGCATTCTGCTGAACCGCTCTCCGATCTTCTTAGATAGGGATAGTGAGmUT
C2	MDR1_86	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGATGTGCTGTATTATCAGmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTCATGAAATGTCCATCTTGAmUA
D2	Pf_Pf3D7_1460900-1_Pf3D7_14	ACACTCTTTCCCTACACGACGCTCTCCGATCTCCCCAAAAGACAATAAGAAAmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGCAAGAGTACTGTTTTATTTmCG
E2	Pf3D7_01_v3_145515_294I_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTTCCGAGTTTTAAGTGAATGmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTAGTTTGTGATGAAGAAAGmGA
F2	Pf3D7_01_v3_180554_D714N	ACACTCTTTCCCTACACGACGCTCTCCGATCTACACCCAGAAATATAAATGGAmAC	TCGGCATTCTGCTGAACCGCTCTCCGATCTTCAAACACGCTTATTACAmAT
G2	Pf3D7_01_v3_535211_2521F_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTAAATAATAAGAACAACGATGmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCACTCATCAATTAATTGTTmAT
H2	Pf3D7_02_v3_470013_G75E_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTAATATGTTGTTGTGTGGTmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGGCATGATAAACCTTTTmAT
A3	Pf3D7_02_v3_714480_D258G	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGAGGTGATCACTATGTTTAmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTCATCCAAGGTTTTGTCTTmAC
B3	Pf3D7_03_v3_155697_150P_B	ACACTCTTTCCCTACACGACGCTCTCCGATCTACGTCTCATTACCAAATTCAmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTTATCAACTGAAGGTACAAmCG
C3	Pf3D7_03_v3_656861_129V_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTATTATTATCATCAATGGTmUG	TCGGCATTCTGCTGAACCGCTCTCCGATCTAGTATGAAGGTGATGATACmUT
D3	Pf3D7_04_v3_1037656_2776I	ACACTCTTTCCCTACACGACGCTCTCCGATCTTAAATGCTGAAGATGAACCAmAG	TCGGCATTCTGCTGAACCGCTCTCCGATCTATGGTAAATAAAGTGAGTGTmGC
E3	Pf3D7_04_v3_139051_K438N_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTAAATTCTAGTAGTAAATACAmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGCATATTCTTTAAATCGTTTmCT
F3	Pf3D7_04_v3_426436_D560A	ACACTCTTTCCCTACACGACGCTCTCCGATCTGGGAGTTGGCTGTAATAATTmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCTCCATGACAAGGAAGAATmUA
G3	Pf3D7_04_v3_531138_A992E_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTTCTCATCCATATGTATCCACmUC	TCGGCATTCTGCTGAACCGCTCTCCGATCTACAAAAGCAGCAATACACATmAT
H3	Pf3D7_04_v3_881571_1081R_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTACAAGGTACATTATTATGGAAmAC	TCGGCATTCTGCTGAACCGCTCTCCGATCTAATATTATATGCAACACCCAmCC
A4	Pf3D7_05_v3_1204155_1338I_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTAGATGAGTGGAAATAAATTTCAmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGTCATATTGGAAGGATCAAAmCT
B4	Pf3D7_05_v3_172801_E218K_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTATTATTGGACGTTCAAGAAAmUG	TCGGCATTCTGCTGAACCGCTCTCCGATCTGGTTCCAATACTAATATGGATmGG
C4	Pf3D7_05_v3_369740_907L_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTGCCTCTATAATCATAGTGAAmUG	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTGTCAATTATTATTGATCTTmUT
D4	Pf3D7_06_v3_1289212_125T_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTGGGAACTATTTAGATTACAmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTTCAAATTTCCCTTACTmGT
E4	Pf3D7_06_v3_900278_P696S_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGTTTCATCTTTATTTCAACAmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGTCATTGGATTTAATTCmUC
F4	Pf3D7_07_v3_1066698_G483S_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTTATTATGTTTAAAGGAAGmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGTTTCAAAATATCTTCATmGT

Well	Name	Forward sequence	Reverse Sequence
G4	Pf3D7_07_v3_1256331_L321F_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTGGTAATATAAATATGTATGATGGAAmGG	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTCACTTTCTTGAACCTACAmCC
H4	Pf3D7_07_v3_1358910	ACACTCTTTCCCTACAGACGCTCTCCGATCTGGAGTTATTCAACTAAGmGC	TCGGCATTCTGCTGAACCGCTCTCCGATCTAGCAAAAATTCTCCATGCMCA
A5	Pf3D7_07_v3_619957_675R_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTCCTATGCTAAAAGTGAACGAmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGTGTATATGGAAGAAAAGmAT
B5	Pf3D7_07_v3_704373_389E	ACACTCTTTCCCTACAGACGCTCTCCGATCTGAAGGATTAAGGAGAAAACAmAG	TCGGCATTCTGCTGAACCGCTCTCCGATCTATCCCTATCTCTTTCTCTTmCC
C5	Pf3D7_08_v3_1056829_L474I_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTGACAATATGGCTAGTAACAmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTATCATCATTGTCTAAAGCTTmCG
D5	Pf3D7_08_v3_1314831_1342K_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTAGAAGATTTTAAAGAAGAAGAmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTCCCTTATATCCATTCTTmUT
E5	Pf3D7_08_v3_150033_1315I_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTATCAACAAGACGTTTCTGATmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTGGTATAACACGTTCCAATAmUT
F5	Pf3D7_08_v3_399774_421K_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTTTCATCATTGTGTGAATGmUC	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTGCTCTTGTGAATGAAAGTmUA
G5	Pf3D7_08_v3_417335_R244K_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTAGGATCATATCTTCGGACTTmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTACAGTAGAGCAAACAAAAGAmAA
H5	Pf3D7_08_v3_549993	ACACTCTTTCCCTACAGACGCTCTCCGATCTAGTGCTTGTACAAAATAATmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTATCCTATAACGAGGTATTmCA
A6	Pf3D7_09_v3_452690_1018I	ACACTCTTTCCCTACAGACGCTCTCCGATCTCCAGGAACCATACTTTTGTmUA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATATTTGACCTGCTTCAAmUG
B6	Pf3D7_09_v3_900277_1534E_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTGATTTGGAATAAAGTATGmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGCCATGTTTAAATGCTTmUG
C6	Pf3D7_10_v3_1383789_N114H	ACACTCTTTCCCTACAGACGCTCTCCGATCTTTAGGTTGGTTAGAATGGAmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTATGAACCAACTTTCTTTACmAT
D6	Pf3D7_10_v3_1386850_927K_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTGGTAAATATATTGAAGACGTmGG	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGAGGACAAGGAAAATAATAmCA
E6	Pf3D7_10_v3_361684	ACACTCTTTCCCTACAGACGCTCTCCGATCTACTTCATCAGCATTTTCAAmCC	TCGGCATTCTGCTGAACCGCTCTCCGATCTCTCATTATTAGGTGGTmGT
F6	Pf3D7_11_v3_1006911_D124E_B	ACACTCTTTCCCTACAGACGCTCTCCGATCTTAAATTTTGC AAAATAGCGTmGG	TCGGCATTCTGCTGAACCGCTCTCCGATCTCGTTTTGGTCATTATGTCTmUA
G6	Pf3D7_11_v3_1020397_G700E_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTTATAATGCATGTGTACCTmUC	TCGGCATTCTGCTGAACCGCTCTCCGATCTACAGAACATTAAACACAACmCA
H6	Pf3D7_11_v3_1295068_E405K	ACACTCTTTCCCTACAGACGCTCTCCGATCTTAATTTGGGAAAATATAAmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTCCCTTCATCATCATTATCATmUT
A7	Pf3D7_11_v3_1815412_E765Q	ACACTCTTTCCCTACAGACGCTCTCCGATCTATGAGTTGTTATATTTCATGTmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTAACATATGCTGCAGATTTTGTmUA
B7	Pf3D7_11_v3_1935031_I139L	ACACTCTTTCCCTACAGACGCTCTCCGATCTGGATGTTCTTTTATGAAATCAmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTACGTCAATACAAGAATTAACmCA
C7	Pf3D7_11_v3_477922_H147Y_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTAACAACCTTAGCATTGATTmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATATTGACTTTTCATCATTmGA
D7	Pf3D7_12_v3_1667593_2381N_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTCGCTGCTGAATATACATATmUG	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTCTTTTCTAATTCCTTTmCA
E7	Pf3D7_12_v3_2171901_V140D_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTGAAGATTCTAAGGAACAAAmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCTGATTCCACAAAAGAAGAAmGA
F7	Pf3D7_12_v3_858501_Q469K	ACACTCTTTCCCTACAGACGCTCTCCGATCTAAATAGAAAATTTGCCACATmGC	TCGGCATTCTGCTGAACCGCTCTCCGATCTAGACCCCTGTTTTAGTAAGAmAT
G7	Pf3D7_12_v3_974663	ACACTCTTTCCCTACAGACGCTCTCCGATCTGAGAAGGAAGACCTTGTmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTAGAATCCTTAAGAGACTGTmUT
H7	Pf3D7_13_v3_1056452_1234D	ACACTCTTTCCCTACAGACGCTCTCCGATCTTTTCATCGCAGGAAAATAmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTCGTCAATATGTGTTATTmCA
A8	Pf3D7_13_v3_1419519	ACACTCTTTCCCTACAGACGCTCTCCGATCTTCTTCTTGACCAACATAmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAGTACAAAACAGmCA
B8	Pf3D7_13_v3_1867630_M4911I	ACACTCTTTCCCTACAGACGCTCTCCGATCTCTCCGATTGGTAGTTATACAmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAAGGAGATGTATCATTmGmAC
C8	Pf3D7_13_v3_2377887_2002S_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTATTAAGAAGAATGTGGAGCmUG	TCGGCATTCTGCTGAACCGCTCTCCGATCTGGAACGGAAGGATATTTAmAT
D8	Pf3D7_13_v3_2573828_I1153M	ACACTCTTTCCCTACAGACGCTCTCCGATCTCCTTCATAACATATAGCAGCmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATGACACCAATCGTTAATCmAC
E8	Pf3D7_13_v3_388365_S1236R	ACACTCTTTCCCTACAGACGCTCTCCGATCTGTCTTTGTCTCTTTCTTmUC	TCGGCATTCTGCTGAACCGCTCTCCGATCTATCCAACCTAAAACAGCTTAmUT

Well	Name	Forward sequence	Reverse Sequence
F8	Pf3D7_14_v3_137622_1179V_A	ACACTCTTCCCTACACGACGCTCTCCGATCTACCATTCATATTTCTAGGTTmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTACGAAATTTATGAAGATGAGmCA
G8	Pf3D7_14_v3_1757603_D1365G	ACACTCTTCCCTACACGACGCTCTCCGATCTGTAAATGACGAAGAAGTAGCmUA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTCTCAATATGTAAATCATTCAmCA
H8	Pf3D7_14_v3_2164225_2830S_B	ACACTCTTCCCTACACGACGCTCTCCGATCTGCTACAAGGAGAAATGACAAmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCATCTGTAGAATCGGTAAGmAA
A9	Pf3D7_14_v3_2733656_557C_A	ACACTCTTCCCTACACGACGCTCTCCGATCTTCAAACCGTTTCATCTACAmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTAAAAGATTCAGCATCGCTAmUC
B9	Pf3D7_14_v3_3126219	ACACTCTTCCCTACACGACGCTCTCCGATCTAACAAGAGAATCAGGGCAmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCATCTCTAGGTAATAACGTACmCA
C9	PlasII_ref	ACACTCTTCCCTACACGACGCTCTCCGATCTTGAAAATGCTCTTTTCACCmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCATCTAAATTCTGCAACATTTmAT
D9	PlasIV_ref	ACACTCTTCCCTACACGACGCTCTCCGATCTCGACAATGCTTTATTACAmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCATTCAATCTCTAAAGAATTTmAT

GRC1 = 68 Primers in total.

## APPENDIX 2 – GRC2 PRIMERS

Well	Name	Forward sequence	Reverse Sequence
A1	CRT_271	ACACTCTTTCCCTACACGACGCTCTCCGATCTTTCCAATTGTTCACTTCTmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATTTTACCTCTACGACTGTmGT
B1	CRT_356	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGTTAGTTGTATACAAGGTCmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTACGTTGTACCATCATAAACAmUT
C1	CRT_97	ACACTCTTTCCCTACACGACGCTCTCCGATCTTTTGGCTAAAAGAAGCTTAAAmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGGTAGGTGGAATAGATTmUC
D1	DHFR_108_164	ACACTCTTTCCCTACACGACGCTCTCCGATCTGTTGTAGTTATGGGAAGAAcAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTATTCTTGATAAACACGGAAmCC
E1	DHPS_613	ACACTCTTTCCCTACACGACGCTCTCCGATCTAATGGAATACCTCGTTATAGmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTCATTTTGTGTTTCATCATmGT
F1	EXO_415	ACACTCTTTCCCTACACGACGCTCTCCGATCTAATTTACCTGAAGACGTTAAmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTATCCCATTGATATCTATACmCT
G1	K13_resistance_2	ACACTCTTTCCCTACACGACGCTCTCCGATCTTATTATCAATACCTCCAACAmAC	TCGGCATTCTGCTGAACCGCTCTCCGATCTATCGTATGAAAGCATGGGTmAG
H1	K13_resistance_4	ACACTCTTTCCCTACACGACGCTCTCCGATCTAATTACTTGAAACATACCATmAC	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATAGGTGGATTTGATGGTmUA
A2	K13_resistance_6	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGGGTATAGTTAACGGATTTmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAAATTGTTGATGCAAATATmUG
B2	MDR1_1034_1042	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGCATTTAGTTCAGATGATmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTCCAAACCAATAGGCAAAACmAT
C2	MDR1_1226_1246	ACACTCTTTCCCTACACGACGCTCTCCGATCTTCTGCAGAAGATTATACTGTmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTACTAACACGTTTAACTCTmCC
D2	MDR1_184	ACACTCTTTCCCTACACGACGCTCTCCGATCTTTTACATATGCCAGTTCCCTmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTGCAACAGTTCTTATCCCATmUA
E2	Pf_PF3D7_1318100_Pf3D7_13	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGATTGAGGACAAATTACATmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCATATAGTTGTAGGGGAGGTmAG
F2	Pf_PF3D7_1447900_Pf3D7_14_v	ACACTCTTTCCCTACACGACGCTCTCCGATCTTATCCATCTCATTGGCTTTmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTGAATAGTTGATAGAGGTACCmGA
G2	Pf3D7_01_v3_179347_311G_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTCACATATCCAGCCCTCmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGCTTGTGTTATTAATCTGTmAA
H2	Pf3D7_01_v3_283144_H664D_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTATGGATTAACCGCTGATATGmUC	TCGGCATTCTGCTGAACCGCTCTCCGATCTATCATCATTTCCTCATCATmCG
A3	Pf3D7_02_v3_376222_K1929E_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTGATGTGATTTCTCTACGAACmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTAGGTATAACTTACGAAATTmUT
B3	Pf3D7_02_v3_529709_F487L	ACACTCTTTCCCTACACGACGCTCTCCGATCTAAAAGACAAGAGTACAAAAGmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTACCAGCTATTCTTGATATGmUG
C3	Pf3D7_02_v3_839620_260L	ACACTCTTTCCCTACACGACGCTCTCCGATCTAGAATGATTGAAAATTGCmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAATCGTGTGCATCCATTATmUT
D3	Pf3D7_03_v3_548178_R2L_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGGGTAAGTTACATTCTTCTmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGTTCATTATGATTAATTATATGmCT
E3	Pf3D7_04_v3_1102392_E808D_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTTAGAGGTGTTGATGTTAATATmGG	TCGGCATTCTGCTGAACCGCTCTCCGATCTATCATCAAATATACATCCCCmAC
F3	Pf3D7_04_v3_110442_G285E_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTATTGAAAATTATGATAGTmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTAGGAAAGACAGCAACAATAAmAA
G3	Pf3D7_04_v3_286542_H586N	ACACTCTTTCCCTACACGACGCTCTCCGATCTTATATAATGCGTCTCTGTCCmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTGAAAAGGAAACAATAAGmAT
H3	Pf3D7_04_v3_529500_1477Y_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTAATAGGAGACAGAATTTGATmUC	TCGGCATTCTGCTGAACCGCTCTCCGATCTTAAGAATTCCTTAGGACAGmUA
A4	Pf3D7_04_v3_648101_51V	ACACTCTTTCCCTACACGACGCTCTCCGATCTAATTTTACATTATACCTCCmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGATACAAATTCATATAAATAmCA
B4	Pf3D7_04_v3_891732_R4468S	ACACTCTTTCCCTACACGACGCTCTCCGATCTCATTCTCTAGTTTAAACCmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGAGCAAATTCCTGAGCATTAmAC
C4	Pf3D7_05_v3_350933	ACACTCTTTCCCTACACGACGCTCTCCGATCTAGCATCATTTTGcmUC	TCGGCATTCTGCTGAACCGCTCTCCGATCTATATATAATCAACAACAmGC
D4	Pf3D7_05_v3_796714_396K_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTAGTTAATGAAAAGAACCCAAmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATGGTGTTTTACGTTTAmGG
E4	Pf3D7_06_v3_1282691_803K_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTCAAATGTACAAGTGAAGGATmUA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGCTTCATTTGTAATAAATAmCT
F4	Pf3D7_06_v3_574938_I2934L_A	ACACTCTTTCCCTACACGACGCTCTCCGATCTTGACGATGATGAAAACATGAmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTACATTTATACATGCCAACAmAT

Well	Name	Forward sequence	Reverse Sequence
G4	Pf3D7_07_v3_1044052_686K_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTGGAACACGCTTATCATTmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTATGACAACCAAATGAAGGAmAT
H4	Pf3D7_07_v3_1213486_S543N_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTATTCTCCTAACTTACGTCAmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTAACAAACACATACATAAAACmCC
A5	Pf3D7_07_v3_1308383_G1945R	ACACTCTTTCCCTACAGACGCTCTCCGATCTCCATTACCTTTACCTTTCCmUC	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTGATTATGATGATGATGACmGA
B5	Pf3D7_07_v3_1359218_K388N_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTAGTAAATTTGAATGGCATAcMGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTTCACGCTTTTCTTAATTGmUA
C5	Pf3D7_07_v3_461139	ACACTCTTTCCCTACAGACGCTCTCCGATCTGTATATATCACAACAATATTTAmGA	TCGGCATTCTGCTGAACCGCTCTCCGATCTCTGAGATTATTTAATTCTTCTTmAC
D5	Pf3D7_07_v3_635985_T598A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTTACTGCGCCTATCATTATTmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAAAGGGATGGAGATTACAAmUG
E5	Pf3D7_08_v3_1313202_799F_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTGGAATAAATGGATTGAGAGAmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATCTTGAATCCATCTTTmCC
F5	Pf3D7_08_v3_339406_1283C_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTATTTCAAATATCCAACCGCmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATAGATGGAAAGCAACCTAmCA
G5	Pf3D7_08_v3_413067_1044V_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTAATGTGTACGAGAAACATCmAC	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATAAGTTTCTTCTAAATAGATAmGT
H5	Pf3D7_09_v3_1379145_R398Q_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTCCATTCAATCTTTTCGTTmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTCCAGATGTGTTATGGTATCmAA
A6	Pf3D7_09_v3_163977_403D_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTTCATTCAAGTGATCCAmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGTAACATAATGTTTTGCAACmAC
B6	Pf3D7_09_v3_599655_E654D_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTATCATGACAAAGAAACATGTmGG	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGTAGATTGTTTCATATCTTmAT
C6	Pf3D7_10_v3_1385894_815P_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTCATTCCACATGTACGTAAGmAA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGAGTGATAAACCTGCTTATmCT
D6	Pf3D7_10_v3_317581_311I_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTTATCTTGTCTCATGTCATGmUG	TCGGCATTCTGCTGAACCGCTCTCCGATCTCAGCTTCTATACCTAATGCTmUG
E6	Pf3D7_10_v3_336274_I1677V_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTTGTATAGTCATGATCATCmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTACCTGGAGATTATGTAATTmCA
F6	Pf3D7_11_v3_1018899_1199L_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTCAAAGGTTACAGAATATTmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGTGGTATAAGTAACATAAGTAmGT
G6	Pf3D7_11_v3_1294107_84A_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTATCCAATAGACCACGAmAG	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAGTCTATATTGAGTTCAGmCT
H6	Pf3D7_11_v3_1802201_450S	ACACTCTTTCCCTACAGACGCTCTCCGATCTGTGTATTTCTTTCCCTGmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAGGAGATAGTTTACCTGGTmUT
A7	Pf3D7_11_v3_1935227_R73S_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTGATCTAAGGATGTAATTTmGT	TCGGCATTCTGCTGAACCGCTCTCCGATCTAGGAATCAAGAACTTTAGCmAT
B7	Pf3D7_11_v3_408668	ACACTCTTTCCCTACAGACGCTCTCCGATCTCACAGGATGATGAATATAATAAmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTGACGGACATTCTTTCAATmGC
C7	Pf3D7_11_v3_828596_K240E_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTGTCTTTTATATCTGTTmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGACATACATAATAAAAGTGAmCA
D7	Pf3D7_12_v3_1663492_1014E_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTCGATGCACCGATTTTCTATmAT	TCGGCATTCTGCTGAACCGCTCTCCGATCTATATTTTGAATCAGCGCTAmUT
E7	Pf3D7_12_v3_1934745_241L_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTACCTGGAATATTACACCTTmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTTGACAAATCATGAAAATCAAmGA
F7	Pf3D7_12_v3_857245_E50G_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTAAGATGTCCTATCAAGACGmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTACGACATACTGCTATTTATGmUT
G7	Pf3D7_13_v3_1233218_N277S_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTGGGTTCTACATTTCCAGATmUT	TCGGCATTCTGCTGAACCGCTCTCCGATCTACCAACTTTGTAAGCTGTAmAA
H7	Pf3D7_13_v3_1466422	ACACTCTTTCCCTACAGACGCTCTCCGATCTTACTTGTAACAGTGCATmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTGTATACGTCGTCATAATTmGG
A8	Pf3D7_13_v3_159086_21R	ACACTCTTTCCCTACAGACGCTCTCCGATCTTCTAACAAACACATTTTATTmCT	TCGGCATTCTGCTGAACCGCTCTCCGATCTCTACTTGGTAGGTATACGTmAT
B8	Pf3D7_13_v3_2161975_D252V	ACACTCTTTCCCTACAGACGCTCTCCGATCTTCCAATCCATGTATATTCATmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTTATGAACGAACTGGACAAmGA
C8	Pf3D7_14_v3_107014_215K	ACACTCTTTCCCTACAGACGCTCTCCGATCTCCTATTCCATACAACATCAAmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAAAGGTATACATGGTATATCmAT
D8	Pf3D7_14_v3_2355751_H1589Q	ACACTCTTTCCCTACAGACGCTCTCCGATCTTATCCAACCCATTTAACGmAG	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAATGTTTCGTTATATATGAmCA
E8	Pf3D7_14_v3_2625887_M238I_A	ACACTCTTTCCCTACAGACGCTCTCCGATCTTCTATTTCGATCTGTACCTTmCA	TCGGCATTCTGCTGAACCGCTCTCCGATCTAAAAGCTTGGTTAGAGATTTmCT

Well	Name	Forward sequence	Reverse Sequence
F8	Pf3D7_14_v3_3046108_417V_A	ACACTCTTTCCCTACACGACGCTCTTCCGATCTCAGGAAATGAGATAATTGACmGT	TCGGCATTCTGCTGAACCGCTCTTCCGATCTTTTCTCATTGGGACGGAmUT
G8	Pf3D7_14_v3_438592_N348T_A	ACACTCTTTCCCTACACGACGCTCTTCCGATCTGAGAAAATGAATCTATGCCTmCA	TCGGCATTCTGCTGAACCGCTCTTCCGATCTATGTGGTCAAGTGATGTATCmAT
H8	PlasI_ref	ACACTCTTTCCCTACACGACGCTCTTCCGATCTCGAACAAGCCGTTTTTACCmUT	TCGGCATTCTGCTGAACCGCTCTTCCGATCTCATCTAAACCTTCAAAGAATTTmGT
A9	PlasIII_ref	ACACTCTTTCCCTACACGACGCTCTTCCGATCTTGAACAAGCCGTTTATTCCmAT	TCGGCATTCTGCTGAACCGCTCTTCCGATCTCACTTGCAGATTCTACGAATTTmAT
B9	W1421_Pf3D7_14_v3_289610_T4	ACACTCTTTCCCTACACGACGCTCTTCCGATCTCTAGGTGACCCATTTATGmAG	TCGGCATTCTGCTGAACCGCTCTTCCGATCTTAGCTTTAGCATCATTCAmCG


GRC2 = 66 Primers in total.

### APPENDIX 3 – SPECIATION PRIMERS


Well	Name	Forward sequence	Reverse Sequence
A1	Spec1	ACACTCTTTCCCTACACGACGCTCTTCCGATCTACCATCCAATTTGATTmGG	TCGGCATTCTGCTGAACCGCTCTTCCGATCTCTTGAGGCAGTTTGTmCC
A2	Spec2	ACACTCTTTCCCTACACGACGCTCTTCCGATCTTAGCGTGTATTGTTGCCTmUG	TCGGCATTCTGCTGAACCGCTCTTCCGATCTCCTTAATGTAGTTTCCTCACAmGC

SPEC = 2 Primers in total.

### APPENDIX 4 – PRIMER POOL:




## APPENDIX 5 – PCR\_1 & PCR\_2




**STORAGE** - Once the PCR is complete overnight at 4°C without detriment, move directly to pooling, or plates may be stored at -20°C for up to one week



**APPENDIX 6 – CLEAN-UP & SIZE SELECTION**



## Appendix 7 – Final size selection & Elution



## Appendix 8 – Primer pooling

